

Product Datasheet

ALDH1A2 peptide (orb234636)

Description	ALDH1A2 Peptide
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	ALDH1A2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb235005 (ALDH1A2 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q63639, Q62148, O94788
Entrez	19378, 8854, 116676
Expiration Date	12 months from date of receipt.