

Product Datasheet

BCL2 (phospho-S87) peptide (orb234675)

Description	BCL2 (phospho-S87) peptide
Reactivity	Human
Tested Applications	Blocking
Target	BCL2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb235044 (BCL2 (phospho-S87) antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P10415
Entrez	596
Expiration Date	12 months from date of receipt.