

Product Datasheet

BCL2 peptide (orb234676)

Description	BCL2 peptide
Reactivity	Human, Mouse, Porcine, Rat, Sheep
Tested Applications	Blocking
Target	BCL2L1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb235045 (BCLX (phospho-S62) antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P53563, Q07817, Q64373
Entrez	12048, 598, 24888
Expiration Date	12 months from date of receipt.