

Product Datasheet

CCR5 (pS336) peptide (orb234696)

Description	The peptide is used to anti-CCR5 (pS336) antibody reactivity
Reactivity	Human, Primate
Tested Applications	Blocking
Target	CCR5
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb235065 (CCR5 (phospho-S336) antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P51682, P51681, O08556
Entrez	12774, 1234, 117029
Expiration Date	12 months from date of receipt.