

Product Datasheet

CNK2 peptide (orb234736)

Description	The peptide is used to anti-CNK2 antibody reactivity
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CNKS2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb235105 (CNK2 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q80YA9, Q8WXI2, Q9Z1T4
Entrez	22866, 59322, 245684
Expiration Date	12 months from date of receipt.