



Product Datasheet

DnaJC3 antibody (orb234825)


Description	Rabbit polyclonal antibody to DnaJC3
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human DnaJC3. The exact sequence is proprietary.
Target	DNAJC3
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q9R0T3 , Q91YW3 , Q13217
Entrez	63880 , 5611 , 100037258
Dilution Range	WB: 1:500-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of DnaJC3 expressi...



Immunohistochemical analysis of DnaJC3 S...



Immunofluorescent analysis of DnaJC3 sta...