



Product Datasheet

ACC alpha (phospho-S80) antibody (orb234990)


Description	Rabbit polyclonal antibody to ACC alpha (phospho-S80)
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S80 of human ACC alpha protein. The exact sequence is proprietary.
Target	ACACA
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q13085, P11497, Q5SWU9
Entrez	60581, 31, 107476
Dilution Range	WB: 1:500-1000, IHC-P: 1:100-200
Expiration Date	12 months from date of receipt.



Western blot analysis of ACC alpha (Phos...



Immunohistochemical analysis of ACC alph...



Direct ELISA antibody dose-response curv...