


Product Datasheet

CaMK2 alpha/delta (phospho-T286) antibody (orb235056)


Description

Rabbit polyclonal antibody to CaMK2


Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Primate, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T286 of human CaMK2 alpha/delta protein. The exact sequence is proprietary.
Target	CAMK2A; CAMK2D
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q6PHZ2 , P11798 , Q9UQM7 , P11275 , P15791 , Q13557
Entrez	815 , 24246 , 12322 , 25400 , 817 , 108058
Dilution Range	WB: 1:500-1000, IHC-P: 1:100-200, IF/ICC: 1:100-500
Expiration Date	12 months from date of receipt.



Western blot analysis of CaMK2 alpha/del...



Immunohistochemical analysis of CaMK2 al...



Immunofluorescent analysis of CaMK2 alph...