

Product Datasheet

CDK5 (phospho-Y15) antibody (orb235084)

Description Rabbit polyclonal antibody to CDK5 (phospho-Y15)

Species/Host Rabbit

Reactivity Human, Mouse, Porcine, Rat, Sheep

Conjugation Unconjugated

Tested Applications IF, IH, WB

Immunogen KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding Y15 of human CDK5 protein. The exact sequence is proprietary.

Target CDK5

Preservatives Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Form/Appearance Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Clonality Polyclonal


Source Rabbit

Uniprot ID [Q00535](#), [Q03114](#), [P49615](#)


Entrez [1020](#), [140908](#), [12568](#)

Dilution Range WB: 1:500-1000, IHC-P: 1:100-200, IF/ICC: 1:100-500


Expiration Date 12 months from date of receipt.



Western blot analysis of CDK5 (Phospho-Y15)



Immunohistochemical analysis of CDK5 (Phospho-Y15)



Immunofluorescent analysis of CDK5 (Phospho-Y15)