


## Product Datasheet

### DNALI1 antibody (orb23624)


|                            |                                                                                      |
|----------------------------|--------------------------------------------------------------------------------------|
| <b>Description</b>         | DNALI1 antibody                                                                      |
| <b>Species/Host</b>        | Rabbit                                                                               |
| <b>Reactivity</b>          | Human, Mouse                                                                         |
| <b>Conjugation</b>         | Unconjugated                                                                         |
| <b>Tested Applications</b> | ELISA, IF, IHC, WB                                                                   |
| <b>Immunogen</b>           | Recombinant Human Axonemal dynein light intermediate polypeptide 1 protein (1-258AA) |
| <b>Target</b>              | DNALI1                                                                               |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol                          |
| <b>Form/Appearance</b>     | Liquid                                                                               |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                        |
| <b>Note</b>                | For research use only                                                                |
| <b>Isotype</b>             | IgG                                                                                  |
| <b>Clonality</b>           | Polyclonal                                                                           |
| <b>Uniprot ID</b>          | <b>O14645</b>                                                                        |
| <b>Dilution Range</b>      | WB:1:500-2000, IHC-P:1:20-200, IF:1:50-200                                           |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                      |



Western blot analysis of Mouse gonadal &...



Immunohistochemical staining of human te...



Immunohistochemical staining of human pr...