



Product Datasheet

Caspase-7 antibody (orb238201)


Description	Caspase-7 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	Recombinant Human Caspase-7 protein (24-198AA)
Target	CASP7
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P55210
Dilution Range	WB:1:1000-1:5000, IHC-P:1:200-1:500, IF:1:50-1:200, IP:1:200-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of A549(Lane 1), H...



Immunoprecipitation analysis of HEK293 W...



Immunofluorescence analysis of MCF-7 cel...