

Product Datasheet

Death domain-containing protein CRADD antibody (orb238226)

Description	Death domain-containing protein CRADD antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Death domain-containing protein CRADD protein (1-199AA)
Target	CRADD
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P78560
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa (Lane 1), K...