

Product Datasheet

Dynein light chain 1, cytoplasmic antibody (HRP) (orb238241)

Description	Dynein light chain 1, cytoplasmic antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Dynein light chain 1, cytoplasmic protein (1-89AA)
Target	DYNLL1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P63167
Expiration Date	12 months from date of receipt.