

## Product Datasheet

# Glyceraldehyde-3-phosphate dehydrogenase antibody (HRP) (orb238249)

|                            |                                                                              |
|----------------------------|------------------------------------------------------------------------------|
| <b>Description</b>         | Glyceraldehyde-3-phosphate dehydrogenase antibody (HRP)                      |
| <b>Species/Host</b>        | Rabbit                                                                       |
| <b>Reactivity</b>          | Human                                                                        |
| <b>Conjugation</b>         | HRP                                                                          |
| <b>Tested Applications</b> | ELISA                                                                        |
| <b>Immunogen</b>           | Recombinant Human Glyceraldehyde-3-phosphate dehydrogenase protein (3-335AA) |
| <b>Target</b>              | Ovomucoid                                                                    |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol               |
| <b>Form/Appearance</b>     | Liquid                                                                       |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                |
| <b>Note</b>                | For research use only                                                        |
| <b>Isotype</b>             | IgG                                                                          |
| <b>Clonality</b>           | Polyclonal                                                                   |
| <b>Uniprot ID</b>          | <b>P04406</b>                                                                |
| <b>Dilution Range</b>      | WB:1:1000-5000, IF:1:50-200,<br>IP:1:200-2000                                |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                              |