



Product Datasheet

Nuclear factor NF-kappa-B p105 subunit 1 antibody (orb238254)


Description	Nuclear factor NF-kappa-B p105 subunit 1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, IF, IHC, WB
Immunogen	Recombinant Human Nuclear factor NF-kappa-B p105 subunit protein (1-961AA)
Target	NFKB1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P19838
Dilution Range	WB:1:500-1:5000, IHC-P:1:100-1:1000, IF:1:200-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2(Lane 1), ...



Immunohistochemical staining of human pl...



Immunoprecipitation analysis of HeLa cel...