

Product Datasheet

Nuclear factor NF-kappa-B p105 subunit 1 antibody (Biotin) (orb238255)

Description	Nuclear factor NF-kappa-B p105 subunit 1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Nuclear factor NF-kappa-B p105 subunit protein (1-961AA)
Target	NFKB1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P19838
Expiration Date	12 months from date of receipt.