


Product Datasheet

Programmed cell death protein 6 antibody (Biotin) (orb238259)

Description	Programmed cell death protein 6 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Programmed cell death protein 6 protein (1-191AA)
Target	PDCD6
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	075340
Expiration Date	12 months from date of receipt.



Western blot analysis of (Lane 1): EC109...