

Product Datasheet

Tumor necrosis factor R superfamily member 6 antibody (Biotin) (orb238286)

Description	Tumor necrosis factor R superfamily member 6 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Tumor necrosis factor receptor superfamily member 6 protein (26-171AA)
Target	FAS
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P25445
Expiration Date	12 months from date of receipt.