

Product Datasheet

Ubiquitin-like-conjugating enzyme ATG3 antibody (orb238289)

Description Ubiquitin-like-conjugating enzyme ATG3 antibody

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications ELISA, IF, IHC, WB

Immunogen Recombinant Human Ubiquitin-like-conjugating enzyme ATG3 protein (1-314AA)

Target GFP Mouse Polyclonal

Preservatives Preservative: 0.03% Proclin 300.
Constituents: 50% Glycerol

Form/Appearance Liquid

Storage Upon receipt, store at -20°C or -80°C.
Avoid repeated freeze.

Note For research use only


Isotype IgG

Clonality Polyclonal


Uniprot ID **Q9NT62**

Dilution Range WB:1:500-1:5000, IHC-P:1:1000-1:2000,
IF:1:200-1:500


Expiration Date 12 months from date of receipt.



Western blot analysis of K562 whole cell...



Immunohistochemical staining of human CO...



Immunohistochemical staining of human ki...