

Product Datasheet

Ubiquitin-like-conjugating enzyme ATG3 antibody (FITC) (orb238291)

Description	Ubiquitin-like-conjugating enzyme ATG3 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Ubiquitin-like-conjugating enzyme ATG3 protein (1-314AA)
Target	Rubisco Mouse Polyclonal (Large Chain)
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NT62
Expiration Date	12 months from date of receipt.