



Product Datasheet

Carbonic anhydrase 12 antibody (orb238317)


Description	Carbonic anhydrase 12 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Carbonic anhydrase 12 protein (25-301AA)
Target	CA12
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O43570
Dilution Range	WB:1:500-1:2000, IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse kidney ti...



Immunohistochemical staining of human ki...



Immunohistochemical staining of human pa...