



Product Datasheet

D-beta-hydroxybutyrate dehydrogenase, mitochondrial antibody (orb238349)

Description	D-beta-hydroxybutyrate dehydrogenase, mitochondrial antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human D-beta-hydroxybutyrate dehydrogenase, mitochondrial protein (47-343AA)
Target	BDH1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q02338
Dilution Range	WB:1:500-1:2000, IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG-2 whole ce...



mitochondrial antibody