

Product Datasheet

Disintegrin and metalloproteinase domain-containing protein 15 antibody (HRP) (orb238360)

Description	Disintegrin and metalloproteinase domain-containing protein 15 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Disintegrin and metalloproteinase domain-containing protein 15 protein (207-452AA)
Target	ADAM15
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q13444
Expiration Date	12 months from date of receipt.