



Product Datasheet

Dual specificity mitogen-activated protein kinase kinase 3 antibody (orb238361)


Description	Dual specificity mitogen-activated protein kinase kinase 3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Dual specificity mitogen-activated protein kinase kinase 3 protein (1-347AA)
Target	MAP2K3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P46734
Dilution Range	WB:1:500-1:2000, IHC-P:1:20-1:200, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of Raw264.7 whole ...



Immunohistochemical staining of human sk...



Immunofluorescence analysis of Hela cell...