



Product Datasheet

Myeloperoxidase antibody (orb238409)

Description	Myeloperoxidase antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Myeloperoxidase protein (165-278AA)
Target	MPO
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P05164
Dilution Range	IHC-P:1:500-1:1000, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of Hela cell...



Immunohistochemical staining of human li...