



Product Datasheet

Protein disulfide-isomerase antibody (orb238429)

Description	Protein disulfide-isomerase antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Protein disulfide-isomerase protein (29-508AA)
Target	P4HB
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P07237
Dilution Range	WB:1:500-1:5000, IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human CO...



Western blot analysis of Lane 1 : Mouse ...