

Product Datasheet

EPO antibody (FITC) (orb23844)

Description	EPO antibody (FITC)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	FITC
Immunogen	Recombinant Macaca fascicularis Erythropoietin protein (28-192AA)
Target	EPO
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P07865
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.