

Product Datasheet

Succinyl-CoA:3-ketoacid coenzyme A transferase 1, mitochondrial antibody (Biotin) (orb238442)

Description	Succinyl-CoA:3-ketoacid coenzyme A transferase 1, mitochondrial antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Succinyl-CoA:3-ketoacid coenzyme A transferase 1, mitochondrial protein (40-520AA)
Target	OXCT1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P55809
Expiration Date	12 months from date of receipt.