



Product Datasheet

Transforming growth factor beta-3 antibody (orb238453)

Description	Transforming growth factor beta-3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human Transforming growth factor beta-3 proprotein protein (301-412AA)
Target	TGFB3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P10600
Dilution Range	IHC-P:1:50-1:500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human kidney tissue.



Immunohistochemical staining of human placenta tissue.