

Product Datasheet

Tumor necrosis factor ligand superfamily member 13B antibody (orb238457)

Description	Tumor necrosis factor ligand superfamily member 13B antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Tumor necrosis factor ligand superfamily member 13B protein (68-285AA)
Target	TNFSF13B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y275
Expiration Date	12 months from date of receipt.