



Product Datasheet

UTP-glucose-1-phosphate uridylyltransferase antibody (orb238465)


Description	UTP-glucose-1-phosphate uridylyltransferase antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human UTP--glucose-1-phosphate uridylyltransferase protein (1-497AA)
Target	UGP2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q16851
Dilution Range	WB:1:1000-1:5000, IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse liver tis...



Immunohistochemical staining of human li...



Immunohistochemical staining of human li...