


Product Datasheet

UTP-glucose-1-phosphate uridylyltransferase antibody (Biotin) (orb238466)

Description	UTP-glucose-1-phosphate uridylyltransferase antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human UTP--glucose-1-phosphate uridylyltransferase protein (1-497AA)
Target	UGP2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q16851
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse liver tis...