



Product Datasheet

Vascular endothelial growth factor C antibody (orb238469)


Description	Vascular endothelial growth factor C antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Vascular endothelial growth factor C protein (112-227AA)
Target	VEGFC
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P49767
Dilution Range	WB:1:1000-1:5000, IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of 293T whole cell...



Immunohistochemical staining of human pl...



Immunohistochemical staining of human sp...