

Product Datasheet

alpha-Galactosidase A antibody (FITC) (orb238503)

Description	alpha-Galactosidase A antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Alpha-galactosidase A protein (32-429AA)
Target	GLA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P06280
Expiration Date	12 months from date of receipt.