



Product Datasheet

alpha-Galactosidase A antibody (HRP) (orb238504)

Description	alpha-Galactosidase A antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Alpha-galactosidase A protein (32-429AA)
Target	GLA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P06280
Expiration Date	12 months from date of receipt.



Immunohistochemistry
of paraffin embedde...



Immunohistochemistry
of paraffin embedde...