

Product Datasheet

Endothelial protein C R antibody (Biotin) (orb238542)

Description	Endothelial protein C R antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Endothelial protein C receptor protein (18-210AA)
Target	PROCR
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UNN8
Expiration Date	12 months from date of receipt.