

Product Datasheet

Glucagon-like peptide 1 R antibody (orb238546)

Description	Glucagon-like peptide 1 R antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human Glucagon-like peptide 1 receptor protein (24-145AA)
Target	GLP1R
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P43220
Dilution Range	WB:1:1000-1:5000, IF:1:50-1:200
Expiration Date	12 months from date of receipt.

