



Product Datasheet

IDO antibody (orb238577)


Description	IDO antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Indoleamine 2,3-dioxygenase 1 protein (1-403AA)
Target	IDO1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P14902
Dilution Range	WB:1:500-1:5000, IHC-P:1:20-1:200, IF:1:50-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa(Lane 1), k...



Immunohistochemical staining of human ce...



Immunofluorescence analysis of HepG2 cel...