

Product Datasheet

IDO antibody (Biotin) (orb238578)

Description	IDO antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Indoleamine 2,3-dioxygenase 1 protein (1-403AA)
Target	IDO1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P14902
Expiration Date	12 months from date of receipt.