

Product Datasheet

Junctional adhesion molecule B antibody (orb238590)

Description	Junctional adhesion molecule B antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Junctional adhesion molecule B protein (29-238AA)
Target	JAM2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P57087
Expiration Date	12 months from date of receipt.



Western blot analysis of using Junctiona...