


Product Datasheet

Protein AMBP antibody (orb238633)

Description	Protein AMBP antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Protein AMBP protein (20-203AA)
Target	AMBP
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P02760
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...