



Product Datasheet

Protein S100-A4 antibody (orb238637)


Description	Protein S100-A4 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Protein S100-A4 protein (2-101AA)
Target	S100A4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P26447
Dilution Range	WB:1:500-1:5000, IHC-P:1:500-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of A549 whole cell...



Immunohistochemical staining of human pl...



Immunohistochemical staining of human ns...