

Product Datasheet

Troponin T, cardiac muscle antibody (FITC) (orb238667)

Description	Troponin T, cardiac muscle antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Troponin T, cardiac muscle protein (138-298AA)
Target	TNNT2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P45379
Expiration Date	12 months from date of receipt.