



Product Datasheet

Mesothelin antibody (orb238718)


Description	Mesothelin antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP
Immunogen	Recombinant Human Mesothelin protein (37-598AA)
Target	MSLN
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q13421
Dilution Range	IHC-P:1:200-1:500, IF:1:50-1:200, IP:1:200-1:2000
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Immunoprecipitation analysis of HepG2 wh...



Immunofluorescence staining of Hela cell...