

Product Datasheet

Tumor necrosis factor-inducible gene 6 protein antibody (Biotin) (orb238739)

Description	Tumor necrosis factor-inducible gene 6 protein antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Tumor necrosis factor-inducible gene 6 protein (18-277AA)
Target	TNFAIP6
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P98066
Expiration Date	12 months from date of receipt.