

## Product Datasheet

### Active regulator of SIRT1 antibody (orb238759)

|                            |                                                                  |
|----------------------------|------------------------------------------------------------------|
| <b>Description</b>         | Active regulator of SIRT1 antibody                               |
| <b>Species/Host</b>        | Rabbit                                                           |
| <b>Reactivity</b>          | Human                                                            |
| <b>Conjugation</b>         | Unconjugated                                                     |
| <b>Tested Applications</b> | ELISA, IHC, WB                                                   |
| <b>Immunogen</b>           | Recombinant Human 40S ribosomal protein S19 protein (1-145AA)    |
| <b>Target</b>              | RPS19                                                            |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol   |
| <b>Form/Appearance</b>     | Liquid                                                           |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C.<br>Avoid repeated freeze. |
| <b>Note</b>                | For research use only                                            |
| <b>Isotype</b>             | IgG                                                              |
| <b>Clonality</b>           | Polyclonal                                                       |
| <b>Uniprot ID</b>          | <b>P39019</b>                                                    |
| <b>Dilution Range</b>      | WB:1:1000-1:5000, IHC-P:1:20-1:200                               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                  |