



Product Datasheet

ERG antibody (orb23876)


Description	ERG antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP
Immunogen	Recombinant Human Transcriptional regulator ERG protein (1-486AA)
Target	ERG
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P11308
Dilution Range	IHC-P:1:20-200, IF:1:50-500, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of Rat heart...



Immunofluorescence analysis of PC-3 cell...



Immunohistochemical staining of human ad...