

Product Datasheet

DTW domain-containing protein 1 antibody (orb238783)

Description DTW domain-containing protein 1 antibody

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications ELISA, IF, IHC

Immunogen Recombinant Human DTW domain-containing protein 1 protein (1-109AA)

Target DTWD1

Preservatives Preservative: 0.03% Proclin 300.
Constituents: 50% Glycerol

Form/Appearance Liquid

Storage Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.

Note For research use only


Isotype IgG

Clonality Polyclonal


Uniprot ID [Q8N5C7](#)

Dilution Range IHC-P:1:20-1:200, IF:1:50-1:200

Expiration Date 12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Immunofluorescent analysis of MCF-7 cell...