

Product Datasheet

Leucine-rich repeat-containing protein 40 antibody (orb238791)

Description	Leucine-rich repeat-containing protein 40 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Leucine-rich repeat-containing protein 40 protein (1-602AA)
Target	LRRC40
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9H9A6
Expiration Date	12 months from date of receipt.



Western blot analysis of 293T whole cell...