

Product Datasheet

MORF4 family-associated protein 1-like 1 antibody (orb238794)

Description	MORF4 family-associated protein 1-like 1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human MORF4 family-associated protein 1-like 1 protein (1-127AA)
Target	MRAP1L1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96HT8
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 whole cel...